Materials are our **DNA!**

pro **3d** ure  
medical

**A U D I O**

2 0 2 1

Materials are our **DNA!**

Inhalt

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Contents

03	Introduction
04	OTOFLOW® approved!
05	Gel
06	SysTherm® product family
08	audioprint® 3D printing materials
16	Silicone materials
19	Accessories
20	HAPTIKAL®
22	Devices (scanner, software, 3D printer, cleaning, post curing)

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
... is the guiding principle of the dynamic, innovative company **pro3dure** medical. The company's philosophy, focused on meeting specific customer requirements, is founded on an innate strategic imperative that has been the driving force in researching and developing innovative solutions for audiology from the very beginning. In this, just as in our name, 3D is the central pivot of all our activities.

At **pro3dure**, a wealth of experience and concentrated technical know-how are combined with innovative spirit to meet customer-specific demands promptly. To this end, we provide you with customized and sophisticated technology packages. The focus is on generative manufacturing processes and their embedding in digital process chains. We develop customized material concepts for your digital workflow. Let our team give you an insight into highly innovative material concepts and at the same time an outlook into future technology packages for the production of earmolds. As an example, we would like to mention the **HPR series** under our new umbrella brand **audioprint®**, which was specially developed for 3D printing of hearing protection earmolds. We have set ourselves the goal of learning from and with you and thus innovatively advancing the topic of manufacturing earmolds in all areas. Accordingly, a trusting cooperation and proximity to our customers and cooperation partners is of outstanding importance to us in our daily work. Meet our team at the **pro3dure** booth at EUHA 2021 and discuss with us the production of the earmold of today and tomorrow. We are looking forward to meeting you!

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The keywords "validation" and "process" are used today in different contexts and with various definitions. In the following, the term process validation is understood as the proof that a process fulfills the demands at the process results. **The EN ISO 9000:2015 defines a validation as:** Validation is a process. It uses objective evidence to confirm that the requirements which define an intended use or application have been met. Whenever all requirements have been met, a validated status is established. Validation can be carried out under realistic use conditions or within a simulated use environment. The definition of process in EN ISO 13485:2016 refers to EN ISO 9000:2015. **A process is understood as:** A process is a set of activities that are interrelated or that interact with one another. Processes use resources to transform inputs into outputs. Processes are interconnected because the output from one process often becomes the input for another process.

**The combination of both of these theoretical definitions lead to 2 practical questions with respect to the production of earpieces:** "Am I building the earpiece correctly?" (Verification) and "Am I building the correct earpiece?" (Process Validation). Triggered by an increasing demand of validated processes on the part of notified bodies **pro3dure** developed the **OTOFLOW®** concept. **OTOFLOW®** offers verified process steps which can be built up to validated production processes. By this holistic approach it is secured that not only the interactions of different materials but also the interactions of single process steps are considered (e. g. cleaning, post curing, tumbling and further surface finishing). Only on the basis of an overall analysis you can answer in the future with Yes to the following question: "Am I building the correct earpiece?" The flexible **OTOFLOW®** concept paves the way to a validated production process.

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**RPR-1** is the lightcuring thixotropic gel for the repair and modification of earpieces. The transparent dosing canulas guarantee optimal control and application of the material in order to make repair and venting works more convenient and efficient.

<b>RPR-1</b>	klar-transparent / clear-transparent	3 x 3,5 g	Art.-Nr. / item-no.: A1400100
<b>RPR-1</b>	beige / beige	3 x 3,5 g	Art.-Nr. / item-no.: A1400400
<b>RPR-1</b>	rot-transparent / red-transparent	3 x 3,5 g	Art.-Nr. / item-no.: A1400200
<b>RPR-1</b>	blau-transparent / blue-transparent	3 x 3,5 g	Art.-Nr. / item-no.: A1400300

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The light curing signature gel **SGNR-1** can be used to further individualize your earpieces. Due to its optimized thixotropy finest details of e. g. laser engraving can be filled without leaking or smearing during the working process. **SGNR-1** is available in red, blue and white. You can choose between a start packing (3 x 3 g in red, blue and white) or standard packing (2 x 3 g per colour).

Startpackung / start package	3 x 3 g	Art.-Nr. / item-no.: A4020130
Standardpackung rot / standard package red	2 x 3 g	Art.-Nr. / item-no.: A4020131
Standardpackung blau / standard package blue	2 x 3 g	Art.-Nr. / item-no.: A4020132
Standardpackung weiß / standard package white	2 x 3 g	Art.-Nr. / item-no.: A4020133
Spritzenkanülen / syringes canulas	10 St. / 10 pc.	Art.-Nr. / item-no.: A4020134

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**SysTherm®**  
4:1 Dispenser  
Art.-Nr. / item-no.: A4020128

**From now on you can work systhernatically!** Our claim to innovation and technology leadership naturally encompasses machine-specific and non-specific products developed perfectly in step with our material systems. Because of pro3dure's holistic approach we offer you everything in order to produce earpieces with thermoeffect efficiently and under the highest quality of demands.

"From now on you work systhernatically" can be taken literally. The SysTherm® system is not just the world's only two part component thermoeffect material for earmolds with no toxic labelling, but also a unique watersoluble cast material for the digital workflow based on additive manufactured molds. SysTherm® perfectly combines two of the world's firsts. In addition to outstanding chemical and physical characteristics, the SysTherm® polymer possesses further advantages such as non toxic labelling and is available in 4:1 cartridges. Now the customer does not only receive an earpiece with excellent wearing comfort, also the laboratory staff will have maximum confidence for highest production efficiency and safety while working. Additionally, the first water soluble cast material SysTherm® Cast guarantees a tailored solution that is perfectly adapted to the SysTherm® polymer. SysTherm® Cast's new chemical formulation ensures that in combination with the SysTherm® polymer no separator is necessary. Now, you can avoid breaking or cracking the mold from the earpiece. Microcracks in the earmold become a thing of the past. SysTherm® Cast can also be used with silicone materials. The SysTherm® collection is completed by the separator SysTherm® SPT-1 which can be used with conventional 3D-printing cast materials. In conclusion, SysTherm® features a variety of technical innovations to make the cast process more efficeint, robust and convenient.

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Fig. 1: Thermoeffect of SysTherm®,  
Shore A hardness as a result of temperature.

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## From now on all under one roof ...

The 3D printing resins for audiology applications, the **GR series** from **pro3dure medical**, is now united under the umbrella brand **audioprint®** to make your selection easier and clearer.

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The generative manufacturing of earpieces on the base of digital data is nowadays a well established process in the HA business and represents an important part of different process chains. For this digital approach **pro3dure medical** offers a wide variety of technology packages and materials. Therefore, the light curing resins of the **audioprint®** series cover the whole spectra of materials for generative technologies. You can choose from resins for stereolithography to materials for image projection systems (DLP) with different irradiation sources. We offer customized material configurations on request. The **audioprint® GR-1** is a laser curing premium resin on the base of (meth)acrylates for the production of hard earmolds and IE shells. The material distinguishes itself by excellent UV irradiation stability in combination with a high elongation at break. **audioprint® GR-1** offers high process robustness and fulfills the highest demands on aesthetics without brittleness. The field of application for this high tech resin ranges from SLA technology to also DLP systems with UV sources (≤ 388 nm and 405 nm irradiation sources as **audioprint® GR-10**) – so you will stay flexible during your daily routine!

# audio print<sup>®</sup> GR-11 new solution the next generation

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The new **audio print<sup>®</sup> GR-11** is the 3D printing resin for the additive manufacturing of next-generation earmolds. The development focused in particular on the mechanical characteristics of the material in order to increase the stability of the 3D printed earmolds and hearing aid shells. The exceptional mechanical stability compared to the market competitors can be seen in Fig. 3 and 4. **audio print<sup>®</sup> GR-11** earmolds have almost 10 times higher total fracture work, with a 7-8 times higher maximum stress intensity factor compared to the state of the art. **audio print<sup>®</sup> GR-11** is available in 6 colors and in Q1 22 in Europe.

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\*\* in Anlehnung an/on the basis of these data result from measurements of a representative sample, which were determined within the scope of our quality assurance.

# audio print<sup>®</sup> GR-12.1 GR-12.2 Gussformen cast forms

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The generative manufacturing resin **audio print<sup>®</sup> GR-12.1** can be used particularly advantageous for the production of cast forms in order to realize a digital way for silicon earpieces. By defined diffraction particles a super smooth surface of the cast forms can be realized. The low viscosity, optimized mechanical properties and reduced inhibition layer helps you in the steps cleaning, post curing of the cast form and demolding the silicone earpiece.

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pro3dure's **audio print<sup>®</sup> GR-12.2** generative resin has been specially developed for the casting process of polyurethane blanks. Due to the high green flex modulus of the material thin cast molds can be realized. The slight translucency of **audio print<sup>®</sup> GR-12.2** gives you the ability to observe the casting process in order to avoid trapped air in the final part. The material can be used in a favourable manner with the **SysTherm<sup>®</sup>** separator **SPT-1**.

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
3D printing resins of the **audio print<sup>®</sup> HPR series** are the ideal choice for additive manufacturing of hearing protection earmolds for use in work areas with the highest requirements. Hearing protection made of **audio print<sup>®</sup> HPR-1** features outstanding chemical-physical properties such as the highest tensile strength in its class (Fig. 5). In addition, **audio print<sup>®</sup> HPR-1** earpieces have a reversible thermo effect. The material is "hard" when inserted into the ear canal and then becomes flexible due to body temperature, thus providing a comfortable and natural wearing comfort. This effect increases the acceptance of wearing hearing protection and thus occupational safety. Hearing protection made of **audio print<sup>®</sup> HPR-1** can be optimally sealed with **audio print<sup>®</sup> lacquer L-2 FLX**. **audio print<sup>®</sup> HPR-1** is available in colorless, white, blue, red and 3 neon signal colors.

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- outstanding properties
- highest tensile strength in its class
- reversible thermo effect
- natural wearing comfort
- increase of work safety

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**audio print<sup>®</sup> L-1** is a universal light-curing 1-component lacquer for surface sealing of hard earmoulds and IE shells. Due to its optimized flowing and adhesion characteristics a homogeneous and high gloss finish is realized. The surface is scratch resistant and can be easily cleaned. The newly developed additive system is matched to the demands of the acoustician and guarantees highest yellowing protection against long-term UV-irradiation. The lacquer is available in a broad variety of colors in order to provide additional solutions for IEM manufacturing.

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The **audio print<sup>®</sup> L-2 FLX** lacquer is a light-curing soft touch coating for flexible and hard earmoulds. Due to its new chemical formulation you can realize surfaces with highest wearing comfort. By chemical and mechanical bonding the lacquer shows optimum adherence in order to generate a long lasting surface sealing. **audio print<sup>®</sup> L-2 FLX** is an excellent solution for the finishing of generative manufactured earmolds.

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
**audioprint<sup>®</sup> INR** generative resins have been developed for special demands like work protection and audio applications. All resins of the **audioprint<sup>®</sup> INR series** are characterized by highest UV-stability, extraordinary mechanical properties (e. g. highest elongation at break values) and a broad spectrum of new colors. Due to that four new signal colors are now available for 3D-printing of hearing protection earmolds and new special colors for in-ear monitors like cognac-transparent and light blue-opaque. Accordingly, the **audioprint<sup>®</sup> INR series** is the perfect addition to the **audioprint<sup>®</sup> GR-series** and offers you the possibility to fulfill your customer's demands more individually. Make 3D-printing more colorful!

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**PR-1** is a light-curing pouring resin for the traditional production of earmoulds, IE and IEM shells. High elongation at break and excellent UV-stability in combination with **pro3dure's** curing device **CD-2** comply with all demands on modern light curing technology. This provides user-friendliness and cost effectiveness. The material is available in a broad variety of colors.

<b>PR-1</b> 250 g Flasche / bottle	blau-transparent / blue-transparent Art.-Nr. / item-no.: <b>A1000504</b>
klar-transparent / clear-transparent Art.-Nr. / item-no.: <b>A1000500</b>	beige-opak / beige-opaque Art.-Nr. / item-no.: <b>A1000550</b>
rötlich-transparent / reddish-transparent Art.-Nr. / item-no.: <b>A1000503</b>	orange-transluzent / orange-translucent Art.-Nr. / item-no.: <b>A1000534</b>
rot-transparent / red-transparent Art.-Nr. / item-no.: <b>A1000502</b>	gelb-transluzent / yellow-translucent Art.-Nr. / item-no.: <b>A1000529</b>

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**IS-1** is a universal impression material of premium quality. By creamy initial consistency and optimized thixotropic behaviour a convenient dispensing is guaranteed. Its excellent detail drawing, scanability and ease of use are further highlights. (A-silicone, medium bodied, time in the ear approx. 2:30 - 3 min., end hardness 34 Shore A)

**IS-1** 8 x S 50  
Doppelkartuschen à 2 x 25 ml A+B /  
double cartridges of 2 x 25 ml A+B  
orange / orange  
Art.-Nr. / item-no.: **A0010100**

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The earmould silicones of the **ES-series** offer everything you expect from a modern high-strength silicone material for the production of flexible earpieces. The **ES-series** is available in the Shore A hardnesses of 25, 40 and 60 and a broad spectrum of different colors. Due to the optimized flow characteristics and low dispensing pressure generative fabricated cast forms can be filled without air traps in an efficient and controllable way. By a newly developed crosslinker matrix an extraordinary combination of high elongation and high tear and tensile strength is achieved. This enables you to remove the earpiece from the cast with a reduced tendency of microcuts. The working time in a pressure pot is at 50 °C, 60 min.

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otosil SL-1 is ideally suited for smoothing out unevennesses, which can be transferred to silicon earmolds in the "digital workflow" by building steps in 3D-printed cast molds. The varnish is characterized by optimized flow and wetting properties in combination with exceptional adhesion to silicone surfaces. In addition, it is free of toluene and thus supports you in the laboratory with regard to occupational safety. Otoplastics coated with otosil SL-1 literally "shine" due to their surface quality. otosil SL-1 is available from Q1 22.

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- optimal smoothing
- optimized flow and wetting properties
- exceptional adhesion on silicone surfaces
- free from toluene

otosil<sup>®</sup> SL-1  
250 g Flasche / 250 g bottle  
Art.-Nr. / item-no.: A0034150

CL-1  
1 kg Flasche / 1 kg bottle  
Art.-Nr. / item-no.: U1900001

## one-component lacquer for sealing silicone earmolds

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New, water soluble cleaning agent for the removal of resin residues on additively manufactured molded parts.

- non explosive
- ideal for ultrasonic systems
- highly efficient (up to 5 times more than Isopropanol)
- ultrafiltrable
- compatible with CLD-1

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The newly developed **HAPTIKAL®** approach of **pro3dure** offers you unique production possibilities for earpieces based on 3D printed molds by adding new dimensions to the cast process. By implementing surface textures into the mold on the one hand biomimetic functionalities can be placed specifically on the earpiece. On the other hand completely new designs can be generated with **HAPTIKAL®**. So if you wish to become a qualified **HAPTIKAL®** partner, contact us.

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For your individual needs different scanner models are available:

- **free Stage:** the system – independent design-object with high operating comfort
- **mono Scan:** simple and precise and cost effective: generating 3D scan data of models of ear impression – the in-store solution for audiology
- **duo Scan:** modern, quiet operation and system-independent for daily use in the acoustician's shop
- **ds Production:** robust, compact and precise makes it particularly suitable for continuous use under production conditions

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
An importat building block for the digital workflow is the ear impression scanner. pro3dure works in close collaboration with our partner smart optics in order to realize flexible solutions for your laboratory or shop. There are different scanner models available.

This year's highlight is the **free Stage**. It was developed for quality-conscious audiologists and is the ideal solution for trend-oriented shops. The system-independent scanner for digitizing ear canal impressions with its slim design is easy to set up at any workstation and turns it into an eye-catcher. The contemporary design of the 3D scanner combines an elegant look with a robust construction. The smooth surfaces as well as the wide opening allow easy access.

The scanner's ergonomic shape, obtained by deliberately not using a cover for the housing, simplifies the workflow and provides the operator with the opportunity to impress the customer with the entire scanning process. The innovative touch sensor enables the user to start the scanning process directly at the scanner, saving you valuable time.

The consistently high quality of materials used in smart optics products is always most important.

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

AirPod templates

# CAD/CAM Software für profis

## CAD/CAM software for professionals

Fully in line with our process orientated company philosophy and in order to complete the digital workflow we offer you the CAD/CAM packages of our partner **Cyflex**.

Good earmolds provide a service inconspicuously. Through their shape and their material, they fulfill multiple requirements – acoustics and comfort, among others. You can rely on the Secret Ear System with regards to quality and the satisfaction of your clients. Secret Ear Designer is the high-performance earmold modeling software for the digital manufacture of ear protection, earmolds and in-ear-monitoring systems. The application facilitates a comfortable, intuitive and time-saving analysis for users and distinguishes itself through a high degree of automation. With only one license, BTE as well as ITE earmolds can be modeled. Digital technology offers you higher quality and thus less and uncomplicated remakes. Thanks to high performance, you can obtain higher output with the same number of employees. In addition, you have the advantage of being able to develop your own products and have a tool at your disposal to create your own construction data.

## ■ All earmold styles

With just one Secret Ear Designer license you can design IdO and HdO earmolds, as well as ear protection and in-ear-monitors, including molds for soft earmolds.

## ■ Intuitive interface

Thanks to easier handling, employees are in a position to model earmolds that they have designed independently and in detail after a few hours of training.

## ■ Time savings

On the basis of easily adjustable product templates the required modeling steps can be reduced to an absolute minimum.

## ■ ERP connectivity

Two-way transfer of parameters between order management system and Secret Ear Designer by means of ERP-Link.

## ■ Compatibility with output media

The application is compatible with all 3D printing technologies available on the market as well as with milling centers for subtractive production.

## ■ Flexible shape design


For maximum independence and creative freedom, construction elements can be autonomously imported and/or designed.

Secret Ear Builder provides your lab with a tool for complete print job generation. The intuitive software makes maximum use of the printer's build area and allows fast grouping and automated alignment of the earmolds to be printed. All necessary algorithms are available for the creation of support structures optimized for earmolds. All relevant basic settings – also those of your 3D printer – are preconfigured.


Different attractive license models are available for this CAD/CAM package. Let us find the best fitting solution for your individual demands.

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The combination of the **pro3dure** materials and the proven **ASIGA MAX™** form the backbone of your digital process chain. The ASIGA MAX™ is the world's most advanced digital 3D printer. It features maximum productivity while requiring a minimal physical footprint. With 62 µm, HD print precision and a 385 nm UV irradiation source, the complete portfolio of **pro3dure's audioprint®** materials can be used. This makes the ASIGA MAX™ 3D printer an all-rounder.

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
The new **Asiga PRO family** combines proven industry leading precision, reliability and large build envelopes which deliver high-speed print capability for the professional environment. You have the option to choose between several lateral resolutions and build envelopes to fit your needs. The new PRO family can be integrated perfectly into your own production process delivering the highest quality and efficiency.

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The **Miicraft Prime and Profession series printers** offer a wide range of build platform sizes, wavelengths, and xy resolutions to efficiently meet a wide variety of audiology applications. In addition, the proven Miicraft beamer technology gives you the confidence to cope with changing demands in the daily laboratory routine. Of course, validated manufacturing processes for these 3D printers are available for all **audioprint®** materials.

- Prime 110**  
Art.-Nr. / item-no.: U2017001
- Prime 150**  
Art.-Nr. / item-no.: U2017002
- Profession 250**  
Art.-Nr. / item-no.: U2017003

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The new "cleaning device" **CLD-1** represents another building block for the optimized and controlled manufacturing process of medical devices by generative fabrication. Unlike conventional ultrasonic baths the **CLD-1** generates a vortex of the cleaning agent. As a result, microcracks in the component surfaces are significantly reduced, thus further increasing the part quality in terms of surface finish, appearance and mechanics of your dental objects. A selection of different cleaning programs assists you in perfectly matching the cleaning program to the selected material class and choice of cleaning agent. The cleaning container can be easily removed and changed in order to simplify the manufacturing processes in your lab. Changing the cleaning liquid can be done within seconds in a convenient way.

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The **RLR-1** is a roller mixer for homogenization of 3D printing materials with pigment or filler content. Highly filled 3D printing materials like the **audioprint® INR-1** require a thorough mixing before use, handshaking is often insufficient. This can result in significant colour deviations and even build failures may occur. The **RLR-1** Mixer is constructed for constant operation in order to keep your **pro3dure** 3D materials ready for use at any time with the same consistency.

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
The high tech post-curing unit **CD-2** is tuned to the latest generation of 3D printers. It features top-level productivity, the most current LED irradiation technology, sophisticated sensors and a wide range of smart features. Through the touch panel you can fully control and monitor all relevant irradiation parameters. Sensors measure the protective gas flow and inform the user about the status of the light source. This avoids inhibition layers on generated parts. As a result, the biocompatibility and the optimal mechanical properties of your 3D printed parts are guaranteed. In a profound and unprecedented form the **CD-2** combines compactness and efficiency on the highest level of technology, along with an excellent price/performance ratio.

the high tech curing device

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highly efficient

#### GERMANY

pro3dure medical GmbH  
Am Burgberg 13  
58642 Iserlohn  
Phone: +49 (0)2374 920050-0  
Fax: +49 (0)2374 920050-50

#### USA

pro3dure medical LLC  
Valley View Business Center  
9825 Valley View Road  
Eden Prairie, MN 55344  
Phone: +1 (952) 426 1928  
Fax: +1 (952) 681 7515

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Exclusion of liability: We cannot be held responsible for false statements in this brochure.  
Delivery of our products is effected according to our general terms of delivery and payment. Technical modifications reserved.